

DESRESON[®] CONCENTRATE

Concentrated disinfectant for a wide range of surfaces and medical equipment

Product information:

DESRESON[®] CONCENTRATE is a concentrated disinfection cleaner for the disinfection and cleaning of surfaces and medical equipment:

- any medical device and equipment as well as floors
- surfaces from glass, porcelain, metal, rubber, wood, aluminium, and acrylic glass
- furniture and fittings, prosthesis, and dental moulds
- DESRESON[®] CONCENTRATE is pleasantly perfumed

Active ingredients:

100 g containing:
 Quaternary ammonium compound 9.90 g

Microbiological activity and acting time:

<i>Spectrum</i>	<i>1 min</i>	<i>5 min</i>	<i>15 min</i>	<i>30 min</i>	<i>60 min</i>	<i>240 min</i>
<i>Organic load</i>	<i>without</i>	<i>with</i>	<i>without</i>	<i>with</i>	<i>without</i>	<i>with</i>
Bacteria ¹	10.0%	5.0%	2.0%	1.5%	1.0%	0.75%
Viruses ²						
HBV / HIV	1.0%	5.0%				
HCV	1.0%					
Adenovirus		5.0%				
Vacciniavirus		5.0%				
SARS	10.0%	1.0%				
Fungus	10.0%	5.0%	2.0%	1.5%	1.0%	0.75%

¹ Surface disinfection according to DGHM

² Virucidal against enveloped viruses according to RKI recommendation Bundesgesundheitsblatt 01-2004

³ Disinfection time according to DGHM

- algicidal
- anti-mould

Wipe treated surfaces with plain water after disinfection.

Packing unit:

5 litre canister
 10 litre canister

Remarks:

Expert reports:

Prof. Dr. R. Schubert, Frankfurt/ M.; Prof. Dr. E. Marth, Graz; Dr. J. Steinmann, Bremen.

Tested on material compatibility with acrylic glass by Röhm/ Darmstadt.

Keep away from children. Lacrimatory. In case of eye contact rinse carefully with water and contact a physician.

Biodegradable.

DGHM listed.

UBA Reg.-Nr.: 22420001

CE 0124