

## *esteeer*<sup>®</sup> SEPT

### Alcoholic disinfection wipes for surfaces and medical equipment

- Product information:** Wipes for the disinfection of alcohol resistant small surfaces, equipment, and furnishings for medical and dental practices as well as hospitals.  
The treated surfaces dry residue-free within a short period of time. Notice the necessary time of application.
- aldehyde-free
  - no quaternary compounds, therefore no respiratory irritation
  - highly compatible and sparing for different material
  - available perfumed as *Classic* and *Lemon* as well as *Neutral* (without perfume)

- Active ingredients:** 100 g containing:
- |            |         |
|------------|---------|
| ethanol    | 40,50 g |
| n-propanol | 9,20 g  |

**Microbiological activity and acting time:**

<i>Spectrum</i>	<i>Organic load</i>		<i>Acting time</i>
	<i>without</i>	<i>with</i>	
Bacteria <sup>1)</sup>	x	x	30 seconds
Tuberculosis <sup>2)</sup>	x	x	1 minute
Viruses <sup>3)</sup>	HCV/ HBV/ HIV	x	x
	Adenovirus	x	x
	Rotavirus	x	x
	Vacciniavirus	x	x
	SARS	x	x
Fungus <sup>1)</sup>	x	x	30 seconds

<sup>1)</sup> Surface disinfection (active substance solution) according to DGHM

<sup>2)</sup> *M. terrae*

<sup>3)</sup> Virucidal against enveloped viruses according to RKI recommendation Bundesgesundheitsblatt 01-2004

- Packing unit:** Refill bag/ box with 120 wipes  
Refill bag/ box with 150 wipes

**Remarks:** Expert reports for the solution:

Prof. Dr. H.-P. Werner, Schwerin; Prof. Dr. R. Schubert, Frankfurt/ M.; Dr. J. Steinmann, Bremen; Prof. Dr. A. Kramer, Greifswald.

- Contains alcohol. Check for material compatibility. Do not use for acrylic glass.
- Labelling according to Gefahrstoffverordnung (Ordinance on Hazardous Substances): R11: Highly flammable. 27 °C flash point according to DIN 51755.
- Do not swallow, avoid eye contact.
- Keep away from children.
- Keep away from fire. Do not smoke.

Solution under the name 'PROPANO AF' listed by the DGHM.

CE 0124